

Validation Report #029784

Summary

Antigen	Mannma Binding Protein/mannan Binding Lectin (MBP/MBL)
Catalog number	ABIN367216
Supplier	Cusabio
Supplier catalog number	<u>csb-e09717h</u>
Lot number	Y25184072
Method validated	Enzyme-linked immunosorbent assay
Laboratory	CGIBD Advanced Analytics Core
Validation number	<u>029784</u>
Positive Control	Human serum - expression is ~50 ng/mL
Negative Control	Goat serum (non-reactive species)
Notes	Target protein was detected in the positive control sample and not in the negative control sample as expected.



Validation Date: 08/11/14

Full Methods

ELISA kit

• Antigen: Mannma Binding Protein/mannan Binding Lectin (MBP/MBL)

Catalog number: ABIN367216

Supplier: Cusabio

Supplier catalog number: csb-e09717h

Lot number: Y25184072

Controls

• Positive control: Human serum (Sigma Aldrich, Cat# H6914-20ML, Lot# SLBK2170V)

Negative control: Goat serum (Sigma Aldrich, Cat# G9023-10ML, Lot# SLBH2670V)

Protocol

- 1. All reagents were brought up to room temperature for 30 minutes prior to use. The 1x Wash Buffer was prepared by adding 20 mL of 25x Wash Buffer Concentrate to 480 mL of distilled/deionized water and mixing thoroughly.
- 2. The vial of Standard was reconstituted with 1 mL of Sample Diluent, mixed, and allowed to sit for 15 minutes with gentle agitation.
- 3. The standard curve was prepared by creating a 2-fold dilution series of seven standards (including the original undiluted vial) using Sample Diluent. Sample Diluent alone served as the 0 pg/mL standard.
- 4. The assay plate was removed from the foil pouch and 100 μ L of each standard and sample were added to the appropriate wells, in triplicate. The plate was covered with the adhesive strip provided and incubated for 2 hours at 37 °C.
- 5. Approximately 10 minutes before the incubation ended, a 1x Biotin-antibody solution was prepared by diluting 60 μ L of 100x Biotin-antibody into 5940 μ L of Biotin-antibody Diluent.
- 6. The liquid from each well was removed.
- 7. 100 µL of 1x Biotin-antibody solution was added to each well, and the plate was covered with a new adhesive strip, and incubated for 1 hour at 37 °C.
- 8. Approximately 10 minutes before the incubation ended, a 1x HRP-avidin solution was prepared by diluting 60 μ L of 100x HRP-avidin into 5940 μ L of HRP-avidin Diluent.
- 9. Each well was aspirated and washed, repeating the process two times for a total of three washes. Each well was washed by filling each well with 1x Wash Buffer and letting it stand for 2 minutes. After the last wash, remaining Wash Buffer was removed and the plate was inverted and blotted against clean, absorbent paper towels.
- 10. 100 μ L of 1x HRP-avidin solution was added to each well, the plate was covered with a new adhesive strip, and incubated for 1 hour at 37 °C.
- 11. The aspiration/wash procedure from Step 9 was repeated for an additional 5 washes.
- 12. 90 µL of TMB Substrate was added to each well. The plate was protected from light and incubated for 15-30 minutes at 37 °C, with periodic checking to prevent overdevelopment.
- 13. 50 µL of Stop Solution was added to each well and mixed thoroughly. The optical density (OD) of each well was measured within 5 minutes using a microplate reader set to 450 nm.
- 14. A standard curve was generated by plotting the OD value for each standard on the y-axis against the concentration on the x-axis. A line of best fit through the points on the graph was used to generate an equation to calculate MBP concentrations of the samples based on their average OD values.

Experimental Notes

The 100x HRP-avidin provided contained less than 60 μ L. Other than that, there were no experimental challenges noted.

Figures




Figure 1: MBP standard curve graph and equation.

ayout							
	1	2	3	4	5	6	
A	STD1	STD1	STD1	SPL1:1	SPL1:1	SPL1:1	Well ID
	0	0	0	50	50	50	Conc/Dil
				Human	Human	Human	Name
				serum	serum	serum	reamo
	STD2	STD2	STD2	SPL1:2	SPL1:2	SPL1:2	Well ID
В	1562.5	1562.5	1562.5	100	100	100	Conc/Dil
ь				Human	Human	Human	Name
				serum	serum	serum	rvanio
С	STD3	STD3	STD3	SPL1:3	SPL1:3	SPL1:3	Well ID
	3125	3125	3125	200	200	200	Conc/Dil
C				Human	Human	Human	Name
	er turnelarin in			serum	serum	serum	Name
	STD4	STD4	STD4	SPL1:4	SPL1:4	SPL1:4	Well ID
D	6250	6250	6250	500	500	500	Conc/Dil
D				Human	Human	Human	Name
				serum	serum	serum	rvame
	STD5	STD5	STD5	SPL1:5	SPL1:5	SPL1:5	Well ID
Е	12500	12500	12500	1000	1000	1000	Conc/Dil
_				Human	Human	Human	Name
				serum	serum	serum	
	STD6	STD6	STD6	SPL1:6	SPL1:6	SPL1:6	Well ID
F	25000	25000	25000	2000	2000	2000	Conc/Dil
г				Human	Human	Human	Name
				serum	serum	serum	realine.
	STD7	STD7	STD7	SPL1:7	SPL1:7	SPL1:7	Well ID
_	50000	50000	50000	5000	5000	5000	Conc/Dil
G				Human	Human	Human	Name
				serum	serum	serum	wame
	STD8	STD8	STD8	SPL2	SPL2	SPL2	Well ID
	1.00E+05	1.00E+05	1.00E+05	1	1/20	10	Conc/Dil
н				Goat	Goat	Goat	No.
				serum	serum	serum	Name

Figure 2: Plate layout. Standard concentrations are in pg/mL; serum dilution values indicate their fold change from the undiluted stock.

RAW DW							
	1	2	3	4	5	6	
Α	0.358	0.356	0.415	2.199	2.152	2.066	RAW DW
В	0.475	0.497	0.58	2.341	2.265	2.413	RAW DW
С	0.582	0.612	0.663	1.16	1.124	1.174	RAW.DW
D	0.805	0.844	0.89	0.751	0.732	0.755	RAW DW
E	1.185	1.142	1.156	0.48	0.494	0.465	RAW DW
F	1.516	1.561	1.626	0.316	0.318	0.3	RAW DW
G	2.121	1.964	2.168	0.185	0.188	0.178	RAW DW
н	2.352	2.215	2.672	0.05	0.053	0.045	RAW DW

Figure 3: Raw OD readings of standards and controls.

Conc							
00110	1	2	3	4	5	6	
Α	<0.000	<0.000	398.863	64208.65	59179.21	51274.97	Conc
В	1145.308	1423.618	2506.451	83724.45	72347.48	97206.87	Conc
С	2533.263	2939.912	3651.504	12447.95	11663.73	12760.92	Conc
D	5786.171	6417.095	7188.974	4945.768	4658.726	5006.757	Conc
E	13010.08	12052.21	12359.37	1208.286	1385.477	1019.8	Conc
F	22167.26	23727.23	26156.35	<0.000	<0.000	<0.000	Conc
G	56153.62	43583.4	60828.14	< 0.000	< 0.000	<0.000	Conc
н	85591.02	66056.98	>105000.0 00	<0.000	<0.000	<0.000	Conc
Cono y Dil							
Conc x Dil	-	-	3		1.		
	1	2	3	4	5	6	2
Α				3210433	2958960	2563748	Conc x Dil
В				8372445	7234748	9720687	Conc x Dil
С				2489591	2332746	2552185	Conc x Dil:
D				2472884	2329363	2503379	Conc x Dil
Е				1208286	1385477	1019800	Cond x Dil
F	,			<0.000	<0.000	<0.000	Conc x Dil
G				<0.000	<0.000	<0.000	Conc x Dil
Н				<0.000	<0.000	<0.000	Conc x Dil

Figure 4: MBP concentrations calculated from standard curve formula. Upper panel = uncorrected for dilution; lower panel = corrected for dilution. On average, 3490315 pg/mL of MBP was detected in the positive control (human serum) and 0 pg/mL of MBP was detected in the negative control (goat serum).